

# JUTURNA-vq PTA Balloon Catheter

## TECHNICAL SPECIFICATIONS

Description	PTA Balloon Catheter   Semi-Compliant OTW
Platinum X-ray Marker	One marker at proximal and distal end of balloon
Balloon Material	Nylon – 0.014", 0.018" and 0.018" Dialysis Fistulae Pebax – 0.035"
Hydrophilic Coating	400 mm distal – 0.014", 0.018" and 0.035"

### PTA Balloon Catheter – 0.014"

Recommended Guidewire	0.014" (0.36 mm)
Recommended Introducer Sheath	4 F - 5 F
Nominal Pressure	1.25 mm to 5 mm - 9 bar
Rated Burst Pressure (RBP)	18 bar
Working Catheter Length	120 cm   150 cm
Low Entry Profile	0.45 mm (0.018")

### PTA Balloon Catheter – 0.018"

Recommended Guidewire	0.018" (0.46 mm)
Recommended Introducer Sheath	4 F - 6 F
Nominal Pressure	2 mm to 7 mm - 9 bar
Rated Burst Pressure (RBP)	18 bar
Working Catheter Length	80 cm   135 cm
Low Entry Profile	0.50 mm (0.020")

### PTA Balloon Catheter – 0.035"

Recommended Guidewire	0.035" (0.89 mm)
Recommended Introducer Sheath	5F, 6F, 9F
Nominal Pressure	3 mm to 12 mm - 7 bar
Rated Burst Pressure (RBP)	• Ø 3 mm - 13 bar • Ø 3.5 mm to 12 mm - 10 bar
Working Catheter Length	80 cm   120 cm   150 cm
Low Entry Profile	0.95 mm (0.037")

### PTA Balloon Catheter – Dialysis Fistulae 0.018"

Recommended Guidewire	0.018" (0.46 mm)
Nominal Pressure	4 mm to 7 mm - 9 bar
Rated Burst Pressure (RBP)	• Ø 4.0 mm & 5.0 mm - 22 bar • Ø 6.0 mm x 20 mm & 7.0 mm x 20 mm - 20 bar • Ø 6.0 mm x 40 mm & 7.0 mm x 40 mm - 18 bar
Working Catheter Length	45 cm   80 cm
Low Entry Profile	0.50 mm (0.020")

Non stock sizes require a minimum order of 5 units and will take up to 8 weeks for delivery

CE 0124



[www.tsunamed.com](http://www.tsunamed.com)

QualiMed Innovative Medizinprodukte GmbH  
Boschstraße 16, 21423 Winsen, Germany

Phone +49 4171 6578 0  
Fax +49 4171 6578 11  
Web [www.qualimed.de](http://www.qualimed.de)

Q<sup>3</sup>Medical GROUP  
DINCESE LIMITED

MANUFACTURED BY



BRO\_JUTURNA-vq REV. B

PERIPHERAL  
VASCULAR



# JUTURNA-vq

PTA Balloon Catheter

0.014" | 0.018" | 0.035" | Dialysis Fistulae 0.018"



Low-Profile

Regulatory Status  
CE Mark Approved

# JUTURNA-vq

PTA Balloon Catheter  
0.014" | 0.018" | 0.035"


The **JUTURNA-vq** PTA Balloon Catheter is a **high performance** balloon catheter for peripheral indication. The device features an **ultra-low profile**, semi-compliant balloon combined with a **low profile tip**.

Our focus is to enhance technology by developing devices with Low Entry Tip Profiles, guided by Ergonomic Pushability for your ease of use, especially in complex vessels, Tactile Trackability for enhanced Manoeuvring.

## KEY FEATURES & BENEFITS

- > **Lowest entry profile in the industry** with the addition of a beveled tip to allow for access to the most difficult lesions
- > **Maximum flexibility** for easy navigation to reach the most distal lesions
- > **Fast inflation and deflation times** to maximize efficiency during treatment and reducing the procedure time compared to other industry participants
- > **Increased pushability** for a significant amount of applied force transition from the hub to the tip
- > **Reduced radiation exposure** for clinician, team, and patient due to decreased time for balloon operation
- > **High rated burst pressure** to optimize working range

## JUTURNA-vq PTA Balloon Catheter



### ORDER INFORMATION

	Balloon Length 120 cm		Catheter Length 120 cm		Catheter Length 150 cm	
	Balloon Length (mm)	Catheter Length 120 cm	Catheter Length 120 cm	Catheter Length 150 cm	Balloon Length (mm)	Catheter Length 150 cm
Balloon Ø 1.25 mm	10	12PVQ012010CPO	15PVQ012010CPO	120	12PVQ012120CPO	15PVQ012120CPO
	15	12PVQ012015CPO	15PVQ012015CPO	140	12PVQ012140CPO	15PVQ012140CPO
	20	12PVQ012020CPO	15PVQ012020CPO	150	12PVQ012150CPO	15PVQ012150CPO
	30	12PVQ012030CPO	15PVQ012030CPO	170	12PVQ012170CPO	15PVQ012170CPO
	40	12PVQ012040CPO	15PVQ012040CPO	200	12PVQ012200CPO	15PVQ012200CPO
	60	12PVQ012060CPO	15PVQ012060CPO	220	12PVQ012220CPO	15PVQ012220CPO
	80	12PVQ012080CPO	15PVQ012080CPO	250	12PVQ012250CPO	15PVQ012250CPO
	100	12PVQ012100CPO	15PVQ012100CPO	280	12PVQ012280CPO	15PVQ012280CPO
Minimum Balloon Profile 0.55 mm						
Balloon Ø 1.5 mm	10	12PVQ015010CPO	15PVQ015010CPO	120	12PVQ015120CPO	15PVQ015120CPO
	15	12PVQ015015CPO	15PVQ015015CPO	140	12PVQ015140CPO	15PVQ015140CPO
	20	12PVQ015020CPO	15PVQ015020CPO	150	12PVQ015150CPO	15PVQ015150CPO
	30	12PVQ015030CPO	15PVQ015030CPO	170	12PVQ015170CPO	15PVQ015170CPO
	40	12PVQ015040CPO	15PVQ015040CPO	200	12PVQ015200CPO	15PVQ015200CPO
	60	12PVQ015060CPO	15PVQ015060CPO	220	12PVQ015220CPO	15PVQ015220CPO
	80	12PVQ015080CPO	15PVQ015080CPO	250	12PVQ015250CPO	15PVQ015250CPO
	100	12PVQ015100CPO	15PVQ015100CPO	280	12PVQ015280CPO	15PVQ015280CPO
Minimum Balloon Profile 0.60 mm						
Balloon Ø 2 mm	10	12PVQ020010CPO	15PVQ020010CPO	120	12PVQ020120CPO	15PVQ020120CPO
	15	12PVQ020015CPO	15PVQ020015CPO	140	12PVQ020140CPO	15PVQ020140CPO
	20	12PVQ020020CPO	15PVQ020020CPO	150	12PVQ020150CPO	15PVQ020150CPO
	30	12PVQ020030CPO	15PVQ020030CPO	170	12PVQ020170CPO	15PVQ020170CPO
	40	12PVQ020040CPO	15PVQ020040CPO	200	12PVQ020200CPO	15PVQ020200CPO
	60	12PVQ020060CPO	15PVQ020060CPO	220	12PVQ020220CPO	15PVQ020220CPO
	80	12PVQ020080CPO	15PVQ020080CPO	250	12PVQ020250CPO	15PVQ020250CPO
	100	12PVQ020100CPO	15PVQ020100CPO	280	12PVQ020280CPO	15PVQ020280CPO
Minimum Balloon Profile 0.65 mm						

### Recommended introducer sheath

4 F (1.35 mm) up to Ø 4.5 mm      5 F (1.67 mm) for Ø 5 mm

# JUTURNA-vq PTA Balloon Catheter

0.014"



## ORDER INFORMATION

Balloon Length (mm)	Catheter Length 120 cm		Catheter Length 150 cm	
	Balloon Length (mm)	Catheter Length 120 cm	Catheter Length 150 cm	Catheter Length 150 cm
Balloon Ø 2.5 mm	10	12PVQ025010CPO	15PVQ025010CPO	
	15	12PVQ025015CPO	15PVQ025015CPO	
	20	12PVQ025020CPO	15PVQ025020CPO	
	30	12PVQ025030CPO	15PVQ025030CPO	
	40	12PVQ025040CPO	15PVQ025040CPO	
	60	12PVQ025060CPO	15PVQ025060CPO	
	80	12PVQ025080CPO	15PVQ025080CPO	
100	12PVQ025100CPO	15PVQ025100CPO		
Minimum Balloon Profile 0.70 mm				
Balloon Ø 3 mm	10	12PVQ030010CPO	15PVQ030010CPO	
	15	12PVQ030015CPO	15PVQ030015CPO	
	20	12PVQ030020CPO	15PVQ030020CPO	
	30	12PVQ030030CPO	15PVQ030030CPO	
	40	12PVQ030040CPO	15PVQ030040CPO	
	60	12PVQ030060CPO	15PVQ030060CPO	
	80	12PVQ030080CPO	15PVQ030080CPO	
100	12PVQ030100CPO	15PVQ030100CPO		
Minimum Balloon Profile 0.79 mm				
Balloon Ø 3.5 mm	10	12PVQ035010CPO	15PVQ035010CPO	
	15	12PVQ035015CPO	15PVQ035015CPO	
	20	12PVQ035020CPO	15PVQ035020CPO	
	30	12PVQ035030CPO	15PVQ035030CPO	
	40	12PVQ035040CPO	15PVQ035040CPO	
	60	12PVQ035060CPO	15PVQ035060CPO	
	80	12PVQ035080CPO	15PVQ035080CPO	
100	12PVQ035100CPO	15PVQ035100CPO		
Minimum Balloon Profile 0.85 mm				
Balloon Ø 4 mm	10	12PVQ040010CPO	15PVQ040010CPO	
	15	12PVQ040015CPO	15PVQ040015CPO	
	20	12PVQ040020CPO	15PVQ040020CPO	
	30	12PVQ040030CPO	15PVQ040030CPO	
	40	12PVQ040040CPO	15PVQ040040CPO	
	60	12PVQ040060CPO	15PVQ040060CPO	
	80	12PVQ040080CPO	15PVQ040080CPO	
100	12PVQ040100CPO	15PVQ040100CPO		
Minimum Balloon Profile 0.95 mm				
Balloon Ø 4.5 mm	10	12PVQ045010CPO	15PVQ045010CPO	
	15	12PVQ045015CPO	15PVQ045015CPO	
	20	12PVQ045020CPO	15PVQ045020CPO	
	30	12PVQ045030CPO	15PVQ045030CPO	
	40	12PVQ045040CPO	15PVQ045040CPO	
	60	12PVQ045060CPO	15PVQ045060CPO	
	80	12PVQ045080CPO	15PVQ045080CPO	
100	12PVQ045100CPO	15PVQ045100CPO		
Minimum Balloon Profile 1.00 mm				
Balloon Ø 5 mm	10	12PVQ050010CPO	15PVQ050010CPO	
	15	12PVQ050015CPO	15PVQ050015CPO	
	20	12PVQ050020CPO	15PVQ050020CPO	
	30	12PVQ050030CPO	15PVQ050030CPO	
	40	12PVQ050040CPO	15PVQ050040CPO	
	60	12PVQ050060CPO	15PVQ050060CPO	
	80	12PVQ050080CPO	15PVQ050080CPO	
100	12PVQ050100CPO	15PVQ050100CPO		
Minimum Balloon Profile 1.05 mm				

### Recommended introducer sheath

4 F (1.35 mm) up to Ø 4.5 mm	5 F (1.67 mm) for Ø 5 mm
------------------------------	--------------------------

# JUTURNA-vq PTA Balloon Catheter

0.018"



Competitor A 0.032" Competitor B 0.022" Competitor C 0.024" Competitor D 0.029" JUTURNA-vq 0.020"

## TIP ENTRANCE DIAMETER compatibility comparison



## ORDER INFORMATION

Balloon Length (mm)	Catheter Length 80 cm		Catheter Length 135 cm	
	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm	Catheter Length 135 cm
Balloon Ø 2 mm	10	08PVQ020010BPO	13PVQ020010BPO	
	15	08PVQ020015BPO	13PVQ020015BPO	
	20	08PVQ020020BPO	13PVQ020020BPO	
	30	08PVQ020030BPO	13PVQ020030BPO	
	40	08PVQ020040BPO	13PVQ020040BPO	
	60	08PVQ020060BPO	13PVQ020060BPO	
	80	08PVQ020080BPO	13PVQ020080BPO	
Minimum Balloon Profile 0.79 mm				
Balloon Ø 2.5 mm	10	08PVQ025010BPO	13PVQ025010BPO	
	15	08PVQ025015BPO	13PVQ025015BPO	
	20	08PVQ025020BPO	13PVQ025020BPO	
	30	08PVQ025030BPO	13PVQ025030BPO	
	40	08PVQ025040BPO	13PVQ025040BPO	
	60	08PVQ025060BPO	13PVQ025060BPO	
	80	08PVQ025080BPO	13PVQ025080BPO	
Minimum Balloon Profile 0.81 mm				
Balloon Ø 3 mm	10	08PVQ030010BPO	13PVQ030010BPO	
	15	08PVQ030015BPO	13PVQ030015BPO	
	20	08PVQ030020BPO	13PVQ030020BPO	
	30	08PVQ030030BPO	13PVQ030030BPO	
	40	08PVQ030040BPO	13PVQ030040BPO	
	60	08PVQ030060BPO	13PVQ030060BPO	
	80	08PVQ030080BPO	13PVQ030080BPO	
Minimum Balloon Profile 0.85 mm				
Balloon Ø 3.5 mm	10	08PVQ035010BPO	13PVQ035010BPO	
	15	08PVQ035015BPO	13PVQ035015BPO	
	20	08PVQ035020BPO	13PVQ035020BPO	
	30	08PVQ035030BPO	13PVQ035030BPO	
	40	08PVQ035040BPO	13PVQ035040BPO	
	60	08PVQ035060BPO	13PVQ035060BPO	
	80	08PVQ035080BPO	13PVQ035080BPO	
Minimum Balloon Profile 0.90 mm				

### Recommended introducer sheath

4 F (1.35 mm) up to Ø 4.5 mm	5 F (1.67 mm) Ø 5 mm to 6 mm	6 F (2.0 mm) Ø 6.5 mm to 7 mm
------------------------------	------------------------------	-------------------------------

# JUTURNA-vq PTA Balloon Catheter

0.018"



## ORDER INFORMATION

Balloon Ø 4 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ040010BPO	13PVQ040010BPO
15	08PVQ040015BPO	13PVQ040015BPO	
20	08PVQ040020BPO	13PVQ040020BPO	
30	08PVQ040030BPO	13PVQ040030BPO	
40	08PVQ040040BPO	13PVQ040040BPO	
60	08PVQ040060BPO	13PVQ040060BPO	
80	08PVQ040080BPO	13PVQ040080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ040100BPO
120	08PVQ040120BPO	13PVQ040120BPO
150	08PVQ040150BPO	13PVQ040150BPO
200	08PVQ040200BPO	13PVQ040200BPO
220	08PVQ040220BPO	13PVQ040220BPO
280	08PVQ040280BPO	13PVQ040280BPO

Minimum Balloon Profile 0.95 mm

Balloon Ø 4.5 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ045010BPO	13PVQ045010BPO
15	08PVQ045015BPO	13PVQ045015BPO	
20	08PVQ045020BPO	13PVQ045020BPO	
30	08PVQ045030BPO	13PVQ045030BPO	
40	08PVQ045040BPO	13PVQ045040BPO	
60	08PVQ045060BPO	13PVQ045060BPO	
80	08PVQ045080BPO	13PVQ045080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ045100BPO
120	08PVQ045120BPO	13PVQ045120BPO
150	08PVQ045150BPO	13PVQ045150BPO
200	08PVQ045200BPO	13PVQ045200BPO
220	08PVQ045220BPO	13PVQ045220BPO
280	08PVQ045280BPO	13PVQ045280BPO

Minimum Balloon Profile 1.05 mm

Balloon Ø 5 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ050010BPO	13PVQ050010BPO
15	08PVQ050015BPO	13PVQ050015BPO	
20	08PVQ050020BPO	13PVQ050020BPO	
30	08PVQ050030BPO	13PVQ050030BPO	
40	08PVQ050040BPO	13PVQ050040BPO	
60	08PVQ050060BPO	13PVQ050060BPO	
80	08PVQ050080BPO	13PVQ050080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ050100BPO
120	08PVQ050120BPO	13PVQ050120BPO
150	08PVQ050150BPO	13PVQ050150BPO
200	08PVQ050200BPO	13PVQ050200BPO
220	08PVQ050220BPO	13PVQ050220BPO
280	08PVQ050280BPO	13PVQ050280BPO

Minimum Balloon Profile 1.20 mm

Balloon Ø 5.5 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ055010BPO	13PVQ055010BPO
15	08PVQ055015BPO	13PVQ055015BPO	
20	08PVQ055020BPO	13PVQ055020BPO	
30	08PVQ055030BPO	13PVQ055030BPO	
40	08PVQ055040BPO	13PVQ055040BPO	
60	08PVQ055060BPO	13PVQ055060BPO	
80	08PVQ055080BPO	13PVQ055080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ055100BPO
120	08PVQ055120BPO	13PVQ055120BPO
150	08PVQ055150BPO	13PVQ055150BPO
200	08PVQ055200BPO	13PVQ055200BPO
220	08PVQ055220BPO	13PVQ055220BPO
280	08PVQ055280BPO	13PVQ055280BPO

Minimum Balloon Profile 1.30 mm

Balloon Ø 6 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ060010BPO	13PVQ060010BPO
15	08PVQ060015BPO	13PVQ060015BPO	
20	08PVQ060020BPO	13PVQ060020BPO	
30	08PVQ060030BPO	13PVQ060030BPO	
40	08PVQ060040BPO	13PVQ060040BPO	
60	08PVQ060060BPO	13PVQ060060BPO	
80	08PVQ060080BPO	13PVQ060080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ060100BPO
120	08PVQ060120BPO	13PVQ060120BPO
150	08PVQ060150BPO	13PVQ060150BPO
200	08PVQ060200BPO	13PVQ060200BPO
220	08PVQ060220BPO	13PVQ060220BPO
280	08PVQ060280BPO	13PVQ060280BPO

Minimum Balloon Profile 1.35 mm

Balloon Ø 6.5 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ065010BPO	13PVQ065010BPO
15	08PVQ065015BPO	13PVQ065015BPO	
20	08PVQ065020BPO	13PVQ065020BPO	
30	08PVQ065030BPO	13PVQ065030BPO	
40	08PVQ065040BPO	13PVQ065040BPO	
60	08PVQ065060BPO	13PVQ065060BPO	
80	08PVQ065080BPO	13PVQ065080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ065100BPO
120	08PVQ065120BPO	13PVQ065120BPO
150	08PVQ065150BPO	13PVQ065150BPO
200	08PVQ065200BPO	13PVQ065200BPO
220	08PVQ065220BPO	13PVQ065220BPO
280	08PVQ065280BPO	13PVQ065280BPO

Minimum Balloon Profile 1.40 mm

Balloon Ø 7 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	10	08PVQ070010BPO	13PVQ070010BPO
15	08PVQ070015BPO	13PVQ070015BPO	
20	08PVQ070020BPO	13PVQ070020BPO	
30	08PVQ070030BPO	13PVQ070030BPO	
40	08PVQ070040BPO	13PVQ070040BPO	
60	08PVQ070060BPO	13PVQ070060BPO	
80	08PVQ070080BPO	13PVQ070080BPO	

Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 135 cm
	100	08PVQ070100BPO
120	08PVQ070120BPO	13PVQ070120BPO
150	08PVQ070150BPO	13PVQ070150BPO
200	08PVQ070200BPO	13PVQ070200BPO
220	08PVQ070220BPO	13PVQ070220BPO
280	08PVQ070280BPO	13PVQ070280BPO

Minimum Balloon Profile 1.40 mm

### Recommended introducer sheath

4 F (1.35 mm)	up to Ø 4.5 mm	5 F (1.67 mm)	Ø 5 mm to 6 mm	6 F (2.0 mm)	Ø 6.5 mm to 7 mm
---------------	----------------	---------------	----------------	--------------	------------------

# JUTURNA-vq PTA Balloon Catheter

0.035"



We provide a diverse range of PTA catheters with **over 1000 sizes** to meet your requirements with balloon lengths up to 280 mm.

Competitor A	Competitor B	Competitor C	JUTURNA-vq
0.048"	0.040"	0.048"	0.037"

**TIP ENTRANCE DIAMETER compatibility comparison**

## ORDER INFORMATION

Balloon Length (mm)	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	10	08PVQ030010APO	12PVQ030010APO	15PVQ030010APO
15	08PVQ030015APO	12PVQ030015APO	15PVQ030015APO	
20	08PVQ030020APO	12PVQ030020APO	15PVQ030020APO	
30	08PVQ030030APO	12PVQ030030APO	15PVQ030030APO	
40	08PVQ030040APO	12PVQ030040APO	15PVQ030040APO	
60	08PVQ030060APO	12PVQ030060APO	15PVQ030060APO	
80	08PVQ030080APO	12PVQ030080APO	15PVQ030080APO	

Balloon Length (mm)	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	100	08PVQ030100APO	12PVQ030100APO	15PVQ030100APO
120	08PVQ030120APO	12PVQ030120APO	15PVQ030120APO	
150	08PVQ030150APO	12PVQ030150APO	15PVQ030150APO	
200	08PVQ030200APO	12PVQ030200APO	15PVQ030200APO	
220	08PVQ030220APO	12PVQ030220APO	15PVQ030220APO	
280	08PVQ030280APO	12PVQ030280APO	15PVQ030280APO	

Minimum Balloon Profile 1.11 mm

Balloon Ø 3.5 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	10	08PVQ035010APO	12PVQ035010APO	15PVQ035010APO
15	08PVQ035015APO	12PVQ035015APO	15PVQ035015APO	
20	08PVQ035020APO	12PVQ035020APO	15PVQ035020APO	
30	08PVQ035030APO	12PVQ035030APO	15PVQ035030APO	
40	08PVQ035040APO	12PVQ035040APO	15PVQ035040APO	
60	08PVQ035060APO	12PVQ035060APO	15PVQ035060APO	
80	08PVQ035080APO	12PVQ035080APO	15PVQ035080APO	

Balloon Length (mm)	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	100	08PVQ035100APO	12PVQ035100APO	15PVQ035100APO
120	08PVQ035120APO	12PVQ035120APO	15PVQ035120APO	
150	08PVQ035150APO	12PVQ035150APO	15PVQ035150APO	
200	08PVQ035200APO	12PVQ035200APO	15PVQ035200APO	
220	08PVQ035220APO	12PVQ035220APO	15PVQ035220APO	
280	08PVQ035280APO	12PVQ035280APO	15PVQ035280APO	

Minimum Balloon Profile 1.12 mm

Balloon Ø 4 mm	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	10	08PVQ040010APO	12PVQ040010APO	15PVQ040010APO
15	08PVQ040015APO	12PVQ040015APO	15PVQ040015APO	
20	08PVQ040020APO	12PVQ040020APO	15PVQ040020APO	
30	08PVQ040030APO	12PVQ040030APO	15PVQ040030APO	
40	08PVQ040040APO	12PVQ040040APO	15PVQ040040APO	
60	08PVQ040060APO	12PVQ040060APO	15PVQ040060APO	
80	08PVQ040080APO	12PVQ040080APO	15PVQ040080APO	

Balloon Length (mm)	Balloon Length (mm)	Catheter Length 80 cm	Catheter Length 120 cm	Catheter Length 150 cm
	100	08PVQ040100APO	12PVQ040100APO	15PVQ040100APO
120	08PVQ040120APO	12PVQ040120APO	15PVQ040120APO	
150	08PVQ040150APO	12PVQ040150APO	15PVQ040150APO	
200	08PVQ040200APO	12PVQ040200APO	15PVQ040200APO	
220	08PVQ040220APO	12PVQ040220APO	15PVQ040220APO	
280	08PVQ040280APO	12PVQ040280APO	15PVQ040280APO	

Minimum Balloon Profile 1.13 mm

### Recommended introducer sheath

5 F (1.67 mm)	up to Ø 5 mm	6 F (2.0 mm)	Ø 5.5 mm to 10.0 mm	9 F (3.0 mm)	for Ø 12 mm
---------------	--------------	--------------	---------------------	--------------	-------------



# JUTURNA-vq

**PTA Balloon Catheter  
Dialysis Fistulae 0.018"**

The **JUTURNA-vq** Dialysis Fistulae Catheter is a **high performance** balloon catheter for peripheral indication. The device features an **ultra-low profile**, semi-compliant balloon combined with a **low profile tip**.


Our focus is to enhance technology by developing devices with Low Entry Tip Profiles, guided by Ergonomic Pushability for your ease of use, especially in complex vessels, Tactile Trackability for enhanced Manoeuvring.

**JUTURNA-vq PTA Balloon Catheter  
Dialysis Fistulae**



## KEY FEATURES & BENEFITS

- > **Very low entry profile** with a beveled tip to cross the most difficult lesions
- > **Maximum flexibility** for easy navigation to reach the most distal lesions
- > **Fast inflation and deflation times** to maximize efficiency during treatment and reducing the procedure time
- > **Increased pushability** for a significant force transition from hub to tip
- > **Reduced radiation exposure** for clinician, team, and patient due to decreased intervention time
- > **High rated burst pressure** to optimize working range



### ORDER INFORMATION

	Balloon Length (mm)	Catheter Length 45 cm	Catheter Length 80 cm	Balloon Length (mm)	Catheter Length 45 cm	Catheter Length 80 cm
Balloon Ø 4 mm	20	04PVQ040020BSD	08PVQ040020BSD	40	04PVQ040040BSD	08PVQ040040BSD
Minimum Balloon Profile 0.95 mm						
Balloon Ø 5 mm	20	04PVQ050020BSD	08PVQ050020BSD	40	04PVQ050040BSD	08PVQ050040BSD
Minimum Balloon Profile 1.20 mm						
Balloon Ø 6 mm	20	04PVQ060020BSD	08PVQ060020BSD	40	04PVQ060040BSD	08PVQ060040BSD
Minimum Balloon Profile 1.35 mm						
Balloon Ø 7 mm	20	04PVQ070020BSD	08PVQ070020BSD	40	04PVQ070040BSD	08PVQ070040BSD
Minimum Balloon Profile 1.40 mm						

### Recommended introducer sheath

4 F (1.35 mm)	up to Ø 4.5 mm	5 F (1.67 mm)	Ø 5 mm to 6 mm	6 F (2.0 mm)	Ø 6.5 mm to 7 mm
---------------	----------------	---------------	----------------	--------------	------------------